**DE**

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**SUPRAMID**

**EN**

**Instructions for use, please read carefully!**

The information given in this package leaflet is updated regularly.  
Please read instructions carefully prior to use.

### **Polyamide 6-6/6**

sterile, nonabsorbable  
**suture material**

#### **DESCRIPTION**

SUPRAMID is a pseudomonofilament nonabsorbable surgical suture material of polyamide 6/6, a polymer of hexamethylenediamine and adipic acid. It is coated with polyamide 6, a polymer of  $\epsilon$ -caprolactam, in order to achieve an even surface. SUPRAMID is available either undyed (white) or dyed black (Carbon black, Colour Index 77266).

SUPRAMID is available in various thicknesses and lengths and also as a needle-thread combination. SUPRAMID satisfies the requirements of the European Pharmacopoeia for sterile nonabsorbable suture material.

#### **INDICATIONS**

SUPRAMID is used for skin closure.

#### **USE**

The suture material should be chosen in accordance with the condition of the patient, the experience of the surgeon, the surgical technique, and the size of the wound.

#### **ACTION**

SUPRAMID is used for coaptation of wound edges. It causes only a minimal inflammatory tissue reaction initially and is then slowly and progressively encapsulated by connective tissue. Although SUPRAMID is a nonabsorbable material, hydrolysis *in vivo* can lead to a progressive loss of tensile strength.

#### **CONTRAINDICATIONS**

As a progressive loss of tensile strength may occur after a prolonged period *in vivo*, SUPRAMID should not be used where tensile strength needs to be maintained over a prolonged period.

## **WARNINGS / PRECAUTIONS / INTERACTIONS**

As with all suture material, care must be taken to ensure that the thread is not damaged during handling. In particular, it must not be kinked or crushed by surgical instruments such as needle holders. When tightening the suture always pull on the thread between the needle and the puncture channel. Do not pull the thread too strongly or over sharp objects. When stretching the thread avoid friction with the rubber glove, as this can damage the thread.

In order not to damage the needle, always grasp it 1/3 – 1/2 of the distance from the end to the point. Do not bend the needle, as this leads to a loss of stability. Because of the risk of infection, the user should take particular care not to incur stab wounds when using surgical needles. Used needles must be disposed of correctly (in order to avoid possible risks of infection).

The risk of wound dehiscence varies with the location of the wound and the suture material used, therefore the user should be familiar with the surgical techniques in which SUPRAMID is to be used. As with all foreign bodies, prolonged contact between SUPRAMID and saline solutions can lead to the formation of calculi (urinary and biliary tracts).

SUPRAMID is used where indicated in accordance with standard surgical suturing and knotting techniques and the experience of the user.

## **ADVERSE REACTIONS**

Inflammatory tissue reactions and transient irritation in the area of the wound may occur initially. The foreignbody effect can exacerbate pre-existing tissue inflammation initially.

## **STERILITY**

SUPRAMID is sterilised by ethylene oxide. Already opened or damaged packs must not be used! Do not resterilise!

## **STORAGE CONDITIONS**

Store at not more than 25 °C protected from moisture and direct heat. Do not use after the expiry date!

## DESCRIPTION OF THE SYMBOLS USED ON THE PACKAGE

**REF** = Reference number

**LOT** = Batch number

 = Use by year - month

 = Consult instructions for use

 = Do not reuse

 = Do not resterilize

 = Do not use if package is damaged

**STERILE EO** = Sterilised using ethylene oxide

 = Upper limit of temperature

 = dyed, monofilament, nonabsorbable

 = dyed, twisted, coated, nonabsorbable

**PA** = Polyamide

 = HIBC code

**CE** 0197 = CE marking and identification number of the notified body. Product conforms to the essential requirements of the Council Directive 93/42/EEC concerning medical devices.

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**RESORBA Medical GmbH**